

Sleep Study- Cloud PAT

Patient Education

Cloud PAT (sleep study) are devices that are used by patients overnight to record breathing and sleep habits. One sticky patch is applied to the top center of the chest. The monitor will record overnight once the patient powers on the device and follows the prompts on the device. The cloud PAT device monitors and records the information and stores it directly on the device. All of the necessary equipment to operate the cloud PAT device for the overnight period, including the sticky patch and watch device, will come in the kit provided. When you are finished with the device please put it in the box it came in and return it to C.A.I. The device will then be downloaded and sent to an off site specialist for reading and interpretation. Your physician will then review those results and provide them to you via your next appointment.

If you have any questions regarding the cloud PAT device or if you are having any technical issues with the monitoring equipment please **contact cloud PAT monitoring first**. If your concern is not resolved in a timely manner or to your satisfaction, please feel free to contact a Cardiovascular Technician at C.A.I. for further assistance.